



T: 866-646-3349
 F: 978-343-9616
 magsales@magmedix.com

26" MRI Non Magnetic Bariatric Wheelchair Parts List

03-WCB26DFB / 03-WCB26DFG / 03-WCB26DLB / 03-WCB26DLG

Left or Right Side is determined while shitting in the Wheelchair

01	Brake Assembly Complete (Left Side)	04-WCBR8006L-01
02	Brake Assembly Complete (Right Side)	04-WCBR8006R-01
03	Rubber Tip for Brake Handle (2 per chair)	04-WCBR8007X-01
04	Push Handle Grip (2 per chair)	04-WCCA8007X-01
05	Kick Tube Cap (2 per chair)	04-WCCA8006X-01
06	Handrim Screws and Spacers, for 24" MRI Wheel, Set of 8 (2 sets per chair)	04-WCHW9213X-01
07	Foot Rest Bumper (2 per chair)	04-WCCA8016X-01
08	Journal Cap for Fork Stem Housing (2 per chair)	04-WCCA3028X-01
09	Foot Plate Tension "O" Ring	04-WCLR7012X-01
10	Foot Plate Screw	04-WCLR7013X-01
11	8" Fork with Nut (2 per chair)	04-WCWH9109X-01
12	Bearing for Fork Stem (4 per chair / 2 each side)	04-WCBE2019X-01
13	5/16" Axle and Nut for 8" Front Wheel (2 per chair)	04-WCWH6108X-01
15	Bearing for 8" Front Wheel (4 per chair / 2 each wheel 24" HD & 26" HD chairs)	04-WCBE1017X-01
17	8" Fork Assembly (24" HD & 26" HD / 2 per chair) (fork, 8" front wheel, bearings, axle and nut)	04-WCWH3108X-01
19	8: Front Wheel for HD Chair with Bearings (2 per chair)	04-WCWH3109X-01
22	24" Rear HD Wheel complete with 5/8" Bearings (2 per chair)	04-WCWH1226X-01
24	5/8" Axle and Nut for 24" Rear Wheel (2 per chair)	04-WCWH6110X-01
27	5/8" Bearing for 24" HD Rear Wheel (4 per chair / 2 each wheel)	04-WCBE1015X-01
30	Hub Cap for 24" HD Rear Wheel (2 Per Chair)	04-WCCA3025X-01
32	Seat Guide (black nylon / 4 per chair)	04-WCCA3026X-01
35	Heel Loop with Hardware for HD Chairs, Left	04-WCAC1518L-01
36	Heel Loop with Hardware for HD Chairs, Right	04-WCAC1518R-01

-CL at the end of the part number means there is a color choice

45	Swing Away Footrest Complete, Left Side (24" & 26" HD chairs)	04-WCFR3107L-01
46	Swing Away Footrest Complete, Right Side 24" & 26" HD chairs)	04-WCFR3107R-01
51	Foot Plate for Foot Rest w/Heel Loop, Left Side (24" & 26" HD chairs)	04-WCFR3112L-01
52	Foot Plate for Foot Rest w/Heel Loop, Right Side (24" & 26" HD chairs)	04-WCFR3112R-01
54	Universal Foot Plate Left or Right (no hole for heal loop, 24" & 25" HD chairs)	04-WCFR3111X-01
55	Calf Pad for Leg Rest Specify Color: 22-Burgundy, 50-Green	04-WCLR31XX-01-CL
62	Removable Left Side Desk Length Arm, (24" & 26" HD chairs) Complete with skirt guard, armrest & button lock Specify Color: 22-Burgundy, 50-Green	04-WCAR2606L-01-CL
63	Removable Right Side Desk Length Arm, (24" & 26" HD chairs) Complete with skirt guard, armrest & button lock Specify Color: 22-Burgundy, 50-Green	04-WCAR2606R-01-CL
65	Desk Length Armrest Pad, (2per chair) Specify Color: 22-Burgundy, 50-Green	04-WCAR2204X-01-CL
66	Skirt Guard, 6" Tall x 13-1/4" Long (2 per chair)	04-WCAR3048X-01
68	Button Lock (2 per chair)	04-WCAR9011X-01
80	Upholstery Back HD, 26" Wide x 18" High Specify Color: 22-Burgundy, 50-Green	04-WCUP1313X-01-CL
81	Upholstery Seat HD, 26" Wide x 20" Deep Specify Color: 22-Burgundy, 50-Green	04-WCUP0510X-01-CL
85	Replacement Safety Belt (24" & 24" HD & 26 HD chairs)	04-WCAC1511X-01
86	Front Foot Rest Holding Pin Repair Kit	04-WCFR3115X-01
91	Pivot Link for 26"W x 16"D, w/ Hardware and Bracket (not shown on drawing)	04-WCCA8022X-01
92	Legrest Complete, Left Side (24" to 26" HD chairs) Specify Color: 22-Burgundy, 50-Green	04-WCLR3107L-01-CL
93	Legrest Complete, Right Side (24" to 26" HD chairs) Specify Color: 22-Burgundy, 50-Green	04-WCLR3107R-01-CL

26" MRI Non Magnetic Bariatric Wheelchair Parts List

03-WCB26DFB / 03-WCB26DFG / 03-WCB26DLB / 03-WCB26DLG

