



T: 866-646-3349  
 F: 978-343-9616  
 magsales@magmedix.com

## 24" MRI Non Magnetic Wheelchair Parts List

03-WCA24DFB / 03-WCA24DFG / 03-WCA24DLB / 03-WCA24DLG

Left or Right Side is determined while shitting in the Wheelchair

01	Brake Assembly Complete (Left Side)	04-WCBR8006L-01
02	Brake Assembly Complete (Right Side)	04-WCBR8006R-01
03	Rubber Tip for Brake Handle (2 per chair)	04-WCBR8007X-01
04	Push Handle Grip (2 per chair)	04-WCCA8007X-01
05	Kick Tube Cap (2 per chair)	04-WCCA8006X-01
06	Handrim Screws and Spacers, for 24" MRI Wheel, Set of 8 (2 sets per chair)	04-WCHW9213X-01
07	Foot Rest Bumper (2 per chair)	04-WCCA8016X-01
08	Journal Cap for Fork Stem Housing (2 per chair)	04-WCCA3028X-01
09	Foot Plate Tension "O" Ring	04-WCLR7012X-01
10	Foot Plate Screw	04-WCLR7013X-01
11	8" Fork with Nut (2 per chair)	04-WCWH9109X-01
12	Bearing for Fork Stem (4 per chair / 2 each side)	04-WCBE2019X-01
13	5/16" Axle and Nut for 8" Front Wheel (2 per chair)	04-WCWH6108X-01
14	Bearing for 8" Front Wheel (4 per chair / 2 each wheel 18" to 24" chairs)	04-WCBE1012X-01
16	8" Fork Assembly (18" to 24" / 2 per chair) (fork, 8" front wheel, bearings, axle and nut)	04-WCWH3106X-01
18	8" Front Wheel with Bearings (2 per chair)	04-WCWH3107X-01
21	24" Rear Wheel complete with 5/8" Bearings (2 per chair)	04-WCWH1225X-01
24	5/8" Axle and Nut for 24" Rear Wheel (2 per chair)	04-WCWH6110X-01
26	5/8" Bearing for 24" Rear Wheel (4 per chair / 2 each wheel)	04-WCBE1013X-01
29	Hub Cap for 24" Rear Wheel (2 per chair)	04-WCCA3024X-01
32	Seat Guide (black nylon / 4 per chair)	04-WCCA3026X-01
34	Heel Loop with Hardware, Left or Right Side	04-WCAC1506X-01

-CL at the end of the part number means there is a color choice

43	Swing Away Footrest Complete, Left Side (22" & 24" chairs)	04-WCFR3105L-01
44	Swing Awa Footrest Complete, Right Side (22" & 24" chairs)	04-WCFR3105R-01
47	Legrest Complete, Left Side (18" to 24" standard chairs)	04-WCLR3104L-01
48	Legrest Complete, Right Side (18" to 24" standard chairs)	04-WCLR3104R-01
49	Foot Plate for Footrest w/ Heel Loop, Left Side (18" to 24" chairs)	04-WCFR3109L-01
50	Foot Plate for Footrest w/ Heel Loop, Right Side (18" to 24" chairs)	04-WCFR3109R-01
53	Universal Foot Plate Left or Right (no hole for heal loop, 18" to 24")	04-WCFR3108X-01
55	Calf Pad for Leg Rest Specify Color: 22-Burgundy, 50-Green	04-WCLR31XX-01-CL
60	Removable Left Side Desk Length Arm, (22" & 24" chairs) Complete with skirt guard, armrest & button lock Specify Color: 22-Burgundy, 50-Green	04-WCAR2607L-01-CL
61	Removable Right Side Desk Length Arm, (22" & 24" chairs) Complete with skirt guard, armrest & button lock Specify Color: 22-Burgundy, 50-Green	04-WCAR2607R-01-CL
65	Desk Length Armrest Pad, (2per chair) Specify Color: 22-Burgundy, 50-Green	04-WCAR2204X-01-CL
66	Skirt Guard, 6" Tall x 13-1/4" Long (2 per chair)	04-WCAR3048X-01
68	Button Lock (2 per chair)	04-WCAR9011X-01
76	Upholstery Back, 24" Wide x 18" High Specify Color: 22-Burgundy, 50-Green	04-WCUP1314X-01-CL
77	Upholstery Seat, 24" Wide x 18" Deep Specify Color: 22-Burgundy, 50-Green	04-WCUP0610X-01-CL
85	Replacement Safety Belt (24" & 24" HD & 26 HD chairs)	04-WCAC1511X-01
86	Front Foot Rest Holding Pin Repair Kit	04-WCFR3115X-01
90	Pivot Link for 24"W x 16"D, w/ Hardware and Bracket (not shown on drawing)	04-WCCA8021X-01

## 24" MRI Non Magnetic Wheelchair Parts List

03-WCA24DFB / 03-WCA24DFG / 03-WCA24DLB / 03-WCA24DLG

